



The new version includes the following upgrades:

- **Completely updated for compatibility with changes to the Genbank website**
 - **This will restore downloading of Genbank sequences and records**
- **Updated for compatibility with Windows Vista**
- **Updated for compatibility with Mac OS 10.5**
- **Numerous speed and usability enhancements**

To install the new version, please read the following directions first.

1. Although we have not had any data losses during testing, to ensure the safety of your existing data, **please back up your user's folder before you start.** This can be found on your hard drive:
C:\Program files\Jellyfish\Users\your_email_username on a PC **OR**
Applications\Jellyfish\Users\your_email_username on a Mac.

Make a copy of the User's folder and save it to your desktop or a flash drive. In the event of a problem, you can place it into the Jellyfish folder after installation to see your files within Jellyfish.

2. **Check which version of Jellyfish you are currently running** (open Jellyfish and look under Help menu at "About Jellyfish"), if you have a version prior to version 3.3 (the current version), we recommend that you uninstall the old version before you install the new one for both Mac and PC users.
3. Go to <http://www.jellyfishsoftware.com> and log in with your email username and password. If you don't know your password, please contact us at jellyfish@fieldscientific.com.
4. Go to the download page and download the appropriate version for your computer – Save/extract the file to your desktop.
5. Double click on the setup file to install new version.
6. Open Jellyfish to check the new version, under Help menu, "About Jellyfish" should say version 3.3.1 and "A Field Scientific Product".